

FootActive

Australian Premium Insoles

Kids



Full-length orthotics for children
with flat feet or fallen arches



Biomechanical support + re-alignment



FootActive Kids Orthotics

Support for little feet!

- full-length orthotic developed by Australian podiatrists
- made of flexible, yet supportive super-EVA
- corrects over-pronation, supports fallen arches and flat feet
- suitable for all children's footwear, especially school and sport shoes
- biomechanical arch support
- deep, stabilising heel cup
- in-built metatarsal support
- shock-absorbing heel and forefoot padding, made of Podflex

FootActive Kids orthotics provide biomechanical support and re-alignment, thereby helping to restore your child's natural walking pattern. By re-aligning the child's feet early on, it will allow the feet to develop more naturally and future biomechanical complaints can be prevented and/or alleviated.

FootActive Kids are useful for children with flat feet, fallen arches or over-pronation (*rolling in of the foot*). The insoles may also help in the relief of **Sever's Disease** (*heel pain*) and growing pains in the knees or legs.

FootActive Kids orthotics are safe to use, however, if your child suffers from a foot deformity, illness or medical condition, we recommend you consult your GP or specialist first, before using FootActive.



Please remove any existing insoles or lining from your child's shoes before inserting FootActive. If required, the orthotics can be trimmed to size with normal scissors, following the trimlines at the toe end. You may clean the orthotics with a warm soapy cloth. Do not soak them or put them in the washing machine or dryer.